



RX Occlusal Contact Definitions

Centric (Occl Contact)	Normal (Tape Pull)	VO1 Visibly Out x 2	VO2 Visibly Out x 2.5	VO3 Visibly Out x 3
Leave specific contact points on occlusal surface as show below.	Tape should pull through. .20mm. <i>Lab tape used is Almore Mylar 12 micron Foil S/S Red</i>	Out .40mm- should appear to be visible but not excessive. (RDL Default)	Out .50mm- should appear as if there is a slight to moderate gap. (Implants)	Out .60mm- should appear excessive.
VO4 Visibly Out x4	VO5 Visibly Out x5			
Out .80mm should appear excessive.	Out 1.00mm should appear extremely excessive.			

